

GENOMICS RESEARCH CORE

NanoString Sample Submission

Researchers are responsible for purchase of codesets and consumables directly from NanoString with shipment directly to Genomics Research Core. However, no shipments can be accepted unless the researcher has first discussed the project with the [Genomics Research Core](#). Please send the package to the attention of the originating Principal Investigator to aid in Genomics Research Core tracking.

Please submit 150 ng at no less than 37 ng/ul for RNA assays or 400 ng at no less than 30 ng/ul for DNA assays to allow for quality control and accurate pipetting. Contact the [Genomics Research Core](#) for advice and assistance if you cannot meet these sample recommendations. Samples should be placed in 96 well plates in row order for efficient processing of batches of 12.

Shipping address:

Debby Hollingshead FAO [Principal Investigator name]
3343 Forbes Avenue, 3rd Floor, Rm 320
Pittsburgh, PA 15260